



Qualatex® Quickie Clips and Ribbon Kit

Seal a 9" or 11" balloon and attach a ribbon in one easy step.

1. Pinch the neck of the inflated balloon closed.




2. Tightly twist the neck several times, then stretch it to ensure a secure seal. Secure one end of the ribbon alongside the balloon neck.



3. Push the twisted balloon neck and ribbon all the way back into the Quickie Clip and snap shut.

Note: Stretching and tightly twisting the balloon neck will ensure a snug fit and the most secure seal.



Contents: 100 Qualatex Quickie Clips and 100 4-foot lengths of White Ribbon